



# CLICKIT ER

EMERGENCY EXTERNAL FIXATION SYSTEM

Universal clamp

Snap and Click Coupling

Sterile

MR Conditional





# Indications

ClickIt ER is an MR Conditional emergency external fixation system for temporary stabilisation of open fractures, unstable fractures and fractures associated with soft tissue injury.

## FEATURES

## ADVANTAGES

Universal clamp

Rod/rod, rod/screw ( $\varnothing$  4.0 -  $\varnothing$  5.0 -  $\varnothing$  6.0 mm) and rod/wire connection

Provisional tightening by hand

Quick application

Snap and Click Coupling

Quick application

Only one instrument required

Simple

Sterile

Ready to use

MR Conditional

Allows up to 3 Tesla magnetic resonance

# Components




## Universal clamp


One clamp for rod/rod, rod/screw, rod/wire connections



Rod/rod connection  $\varnothing$  12 mm



Rod/screw connection  $\varnothing$  12 mm (rod) and  $\varnothing$  4, 5, 6 mm (screw)



Rod/wire connection  $\varnothing$  12 mm (rod) and from 1.8 to 2.5 mm (wire)




### Rods

Diameter 12 mm radiolucent




### Multiscrew clamp

Easy positioning of 2 or 3 screws



### Screws and Wires

Available in different diameter, material and thread lengths




### ER hybrid connecting Post

Allow for hybrid fixator configurations




### “U” rods

Optimal accessory for ankle or tibial pilon




### Arch 5/8

For hybrid configurations



### Half rings

For hybrid configurations



### Radiolucent carbon ring

For hybrid configurations

# Custom Kit ER




Code	Description	Q.ty
5006508ST	Universal clamp	max 6



Code	Description	Q.ty
5006504ST	Multiscrew clamp	max 2



Code	Description
5006510ST	Rod Ø 12 mm L 100 mm
5006512ST	Rod Ø 12 mm L 150 mm
5006514ST	Rod Ø 12 mm L 200 mm
5006516ST	Rod Ø 12 mm L 250 mm
5006518ST	Rod Ø 12 mm L 300 mm
5006520ST	Rod Ø 12 mm L 350 mm




Code	Description
5000592ST	Piercing screw L 250 mm
VA4.150.40	Screw Ø 4.0 mm L 150x40 mm
VA5.150.40	Screw Ø 5.0 mm L 150x40 mm
VA6.250.40	Screw Ø 6.0 mm L 250x40 mm

# Configurations


## UPPER LIMB






## PELVIS



## LOWER LIMB










Ankle Spanning ER Kit (Code 5006600ST)








# Kit ClickIt ER


## ANKLE SPANNING KIT



### 5006600ST

Kit consisting of	Q.ty
Multiscrew clamp	1
Universal clamp	4
Rod Ø 12 mm L 200 mm	1
Rod Ø 12 mm L 300 mm	2
Piercing screw L 250 mm	1
Screw Ø 5.0 mm L 150x40 mm	2


## ANKLE SPANNING KIT PLUS



### 5006600STB

Kit consisting of	Q.ty
Multiscrew clamp	1
Universal clamp	6
Rod Ø 12 mm L 200 mm	2
Rod Ø 12 mm L 350 mm	2
Piercing screw L 250 mm	1
Screw Ø 4.0 mm L 150x40 mm	1
Screw Ø 5.0 mm L 150x40 mm	2


## ELBOW TRAUMA KIT



### 5006650ST01

Kit consisting of	Q.ty
Universal clamp	5
Rod Ø 12 mm L 250 mm	2
Screw Ø 4.0 mm L 150x40 mm	2
Screw Ø 5.0 mm L 150x40 mm	2


## KNEE BRIDGING KIT



### 5006630ST01

Kit consisting of	Q.ty
Universal clamp	6
Rod Ø 12 mm L 250 mm	1
Rod Ø 12 mm L 350 mm	2
Screw Ø 5.0 mm L 150x40 mm	2
Screw Ø 6.0 mm L 250x40 mm	2


## ANKLE SPANNING KIT L



### 5006600STA

Kit consisting of	Q.ty
Multiscrew clamp	1
Universal clamp	4
Rod Ø 12 mm L 200 mm	1
Rod Ø 12 mm L 350 mm	2
Piercing screw L 250 mm	1
Screw Ø 5.0 mm L 150x40 mm	2


## TIBIA KIT



### 5006620ST01

Kit consisting of	Q.ty
Universal clamp	6
Rod Ø 12 mm L 150 mm	2
Rod Ø 12 mm L 200 mm	1
Screw Ø 5.0 mm L 150x40 mm	4


## FEMUR KIT



### 5006640ST01

Kit consisting of	Q.ty
Universal clamp	6
Rod Ø 12 mm L 150 mm	1
Rod Ø 12 mm L 200 mm	2
Screw Ø 6.0 mm L 250x40 mm	4

## HUMERUS KIT



### 5006660ST01

Kit consisting of	Q.ty
Universal clamp	6
Rod Ø 12 mm L 150 mm	3
Screw Ø 5.0 mm L 150x40 mm	4

## ClickIt ER instrument set



## ClickIt ER single use instrument set



### 5006010ST01

Code	Description	Q.ty
5000603M	T-wrench single use	1
5006010M	Cannulated guide Ø 6.0 L 60 - Single use	1
5006000M	Hexagonal L-wrench 3.0 mm - Single use	1
5002074M	Hexagonal L-wrench 5.0 mm - Single use	1

### 5006012ST01

Code	Description	Q.ty
5000603M	T-wrench single use	1
5006012M	Cannulated guide Ø 6.0 L 100 - Single use	1
5006000M	Hexagonal L-wrench 3.0 mm - Single use	1
5002074M	Hexagonal L-wrench 5.0 mm - Single use	1

# CLICKIT ER

## COMPONENTS and ACCESSORIES

Code	Description
5006508ST	Universal clamp
5006504ST	ER multiscrew clamp
5006506ST	ER Hybrid Connecting Post
5006510ST	Rod Ø 12 L 100 mm
5006512ST	Rod Ø 12 L 150 mm
5006514ST	Rod Ø 12 L 200 mm
5006516ST	Rod Ø 12 L 250 mm
5006518ST	Rod Ø 12 L 300 mm
5006520ST	Rod Ø 12 L 350 mm
5006522ST	Rod Ø 12 L 400 mm
5006526ST	Rod Ø 12 L 500 mm
5006528ST	Rod Ø 12 L 600 mm
5.2201A.001	ClickIt ER "U" bar Ø 150
5.2201A.002	ClickIt ER "U" bar Ø 180
5.2201A.003	ClickIt ER "U" bar Ø 210
5006102ST	120 mm Ø 5/8 ring
5006104ST	140 mm Ø 5/8 ring
5006106ST	160 mm Ø 5/8 ring
5006108ST	180 mm Ø 5/8 ring
5006110ST	200 mm Ø 5/8 ring
5006112ST	220 mm Ø 5/8 ring

## CLICKIT ER KIT

5006600ST	Ankle Spanning Kit
5006600STA	Ankle Spanning Kit L
5006600STB	Ankle Spanning Kit Plus
5006610ST01	Ankle Delta Kit
5006620ST01	Tibia Kit
5006630ST01	Knee Bridging Kit
5006640ST01	Femur Kit
5006650ST01	Elbow Trauma Kit
5006650STA	Elbow Trauma Kit Plus
5006660ST01	Humerus Kit

## INSTRUMENTS \*

Code	Description	Q.ty
5006562	Complete ClickIt ER instrument trax	
Comprising:		
5000124	3.0 mm hex screwdriver	1
5000600	T-wrench for screws Ø 4.0 mm	1
5000601	T-wrench for screws Ø 5.0 mm	1
5000602	T-wrench for screws Ø 6.0 mm	1
5000604	T-handle with AO attachment	1
5000605	Hexagonal rod Ø 5.0 mm AO connection	1
5006000	3.0 mm L-shaped Allen key	1
5002074	5.0 mm L-shaped Allen key	1
5000222	Twist drill Ø 3.5 mm	1
5000224	Twist drill Ø 4.5 mm	1
5006010	Drill guide Ø 3.5 mm L 110 mm	1
5006.04	Ø 4.0 mm adapter t-wrench	1
5006.05	Ø 5.0 mm adapter t-wrench	1
5006.06	Ø 6.0 mm adapter t-wrench	1
5006564	ER rod template	1
5006560ER	ClickIt ER sterilization trax	1

\* Implantable and sterile devices are covered by the CE-1936 marking, while non-sterile reusable instruments are covered by the CE marking

## Sterilization

The components of the ClickIt ER system are supplied sterile.

We indicate sterilizing the instruments in saturated steam autoclave according to the following parameters:

**TEMPERATURE:** 134°

**EXPOSURE:** 5 min





**HEADQUARTERS**

Via P. Gobetti, 56r - 16145 Genova  
Tel. +39 010 30801 - Fax +39 010 3080210  
[www.mikai.it](http://www.mikai.it) - [servizio.clienti@mikai.it](mailto:servizio.clienti@mikai.it)

**MANUFACTURING FACILITY**

Via Canestrello, 2 - 36050 Monteviale (VI)  
Tel. +39 0444 950100 - Fax +39 0444 950133  
[www.mikai.it](http://www.mikai.it) - [servizio.clienti@mikai.it](mailto:servizio.clienti@mikai.it)